

MRI ELBOW – SR SPORTS

- Patient in prone position with arm overhead - if not tolerated, put patient in supine position with arm at side
- Axials include volume from above epicondyles to below radial tuberosity
- Coronals include entire elbow, skin to skin, parallel to intercondylar line
- Sagittals include entire elbow, skin to skin, 90° to intercondylar line

SEQUENCE	FOV	Slice	TE	Matrix (min)
Ax T1	14 cm	4 mm	20-40	256x256
Ax PD FS	14 cm	4 mm	20-40	256x256
Cor T1	14 cm	4 mm	min	256x256
Cor PD FS	14 cm	4 mm	20-40	256x256
Cor GRE	14 cm	4 mm	min	256x256
Sag PD FS	14 cm	4 mm	20-40	256x256